

Volume	0,5 mL	1,0 mL	2,0 mL	4,0 mL	5,0 mL
Catalogue no.	MLC05 MLC05N	MLC1 MLC1N	MLC2 MLC2N	MLC4 MLC4N	MLC5 MLC5N
Maximum Volume: Max indicator -20°	0,5	0,9	1,7	3,5	4,5
Max indicator -80°	0,4	0,9	1,5	3,5	4,5
Graduation	No	Yes	Yes	Yes	Yes
Material	Polypropylene (tube and cap)				
Certified RNase, DNase,DNA and Pyrogen free	Yes	Yes	Yes	Yes	Yes
Features	Lip seal closure, Half turn quick seal, with writing area, anti-rotation lid and external thread. Star foot anti-rotation bottom, Self-standing style with conical bottom shape. Withstand temperatures to -196°C. CE Marked for storage of diagnostic samples in the EU. Compliant with IATA for shipping.				
Production Standards & Procedures	Production is GMP oriented. Fulfills the criteria of medical devices according to European Directive 42/93/ECC. MAXLINE Cryo tubes are designed and manufactured in Germany in a clean room environment. Hygienic production is ensured by a sealed area. All produced of raw materials, compliant to applicable regulations and standards.				
Packaging	5x100 pcs. per case / sterile / 10x100 pcs. per case / sterile 5x100 pcs. per case / non-sterile / 10x100 pcs. per case / non-sterile				

Technical Drawing

