

| Volume                                       | 0,5 mL                       | 1,0 mL  | 2,0 mL  | 4,0 mL  | 5,0 mL  |
|--|------------------------------|---------|---------|---------|---------|
| Catalogue no.                                | MLC05-NS                     | MLC1-NS | MLC2-NS | MLC4-NS | MLC5-NS |
| Graduation                                   | No                           | No      | No      | No      | No      |
| Material                                     | Polypropylene (tube and cap) |         |         |         |         |
| Certified RNase, DNase, DNA and Pyrogen free | Yes                          | Yes     | Yes     | Yes     | Yes     |

**Features** Lip seal closure, Half turn quick seal, anti-rotation lid and external thread. Star foot anti-rotation bottom, Self-standing style with conical bottom shape. Withstand temperatures to -196°C. CE Marked for storage of diagnostic samples in the EU. Compliant with IATA for shipping.

**Production Standards & Procedures** Production is GMP oriented. Fulfills the criteria of medical devices according to European Directive 42/93/ECC. MAXLINE Cryo tubes are designed and manufactured in Germany in a clean room environment. Hygienic production is ensured by a sealed area. All produced of raw materials, compliant to applicable regulations and standards.

**Packaging** 2x500 pcs. (tubes) and 1x1000 pcs. (caps) per case / non-sterile  
4x250 pcs. (tubes) and 1x1000 pcs. (caps) per case / non-sterile

**Technical Drawing**

