



Sartoclear Dynamics® Lab

2 steps to Choosing your Cell Culture Harvesting Kit

Simplifying Progress


SARTORIUS


Step 1: Choose Your Volume



Step 1: Choose Your Volume

Step 2: Choose Filtration





Easy to Filter (20 g DE/L)


Your Kit: SDLV-1000-20E0-E



1. 1 L vacuum filtration unit
Sartolab® RF 1,000, 0.22 μ m, PES
2. 2 x 10 g Filter Aid pouch



Step 1: Choose Your Volume



Step 2: Choose Filtration

Hard to Filter (40 g DE/L)

Your Kit: SDLV-1000-40E0-E


X


1. 1 L vacuum filtration unit
Sartolab® RF 1,000, 0.22 μm , PES
2. 4 x 10 g Filter Aid pouch



Step 1: Choose Your Volume

Step 2: Choose Filtration






Easy to Filter (20 g DE/L)

Your Kit: SDLV-0500-10E0-2



1. 500 mL vacuum filtration unit
Sartolab® RF 500, 0.22 µm, PES
2. 10 g Filter Aid pouch per unit





Hard to Filter (40 g DE/L)

Your Kit: SDLV-0500-20E0-E





1. 500 mL vacuum filtration unit Sartolab® RF 500, 0.22 µm, PES
2. 2 x 10 g Filter Aid pouch



Step 1: Choose Your Volume

Step 2: Choose Filtration






Easy to Filter (20 g DE/L)

Your Kit: SDLV-0250-05E0-2



1. 250 mL vacuum filtration unit Sartolab® RF 250, 0.22 μm , PES
2. 5 g Filter Aid pouch per unit





Hard to Filter (40 g DE/L)

Your Kit: SDLV-0250-10E0-2





1. 250 mL vacuum filtration unit Sartolab® RF 250, 0.22 μm , PES
2. 10 g Filter Aid pouch per unit



Step 1: Choose Your Volume

Step 2: Choose Filtration






Easy to Filter (20 g DE/L)

Your Kit: SDLV-0150-02E0-E



1. 150 mL vacuum filtration unit
Sartolab® RF 150, 0.22 µm, PES
2. 2 x 1 g Filter Aid pouch per unit





Hard to Filter (40 g DE/L)

Your Kit: SDLV-0150-05E0-2




1. 150 mL vacuum filtration unit Sartolab® RF 150, 0.22 μm , PES
2. 5 g Filter Aid pouch per unit




Step 1: Choose Your Volume

Step 2: Choose Filtration



Step 1: Choose Your Volume



Step 2: Choose Filtration


Easy to Filter (20 g DE/L)

Your Kit: **SDLV-0150-02E0-E**

1. 150 mL vacuum filtration unit Sartolab® RF 150, 0.22 µm, PES
2. 2 x 1 g Filter Aid pouch per unit



Step 1: Choose Your Volume



Step 2: Choose Filtration

Hard to Filter (40 g DE/L)

Your Kit: SDLV-0150-05E0-2

1. 150 L vacuum filtration unit Sartolab® RF 150, 0.22 μm , PES
2. 5 g Filter Aid pouch per unit




+




X

Step 1: Choose Your Volume

Step 2: Choose Filtration



Step 1: Choose Your Volume



Step 2: Choose Filtration

Easy to Filter (20 g DE/L)

Your Kit: **SDLV-0050-01E0-2**

1. 50 mL vacuum filtration unit
Sartolab® RF 50, 0.22 μ m, PES
2. 1 g Filter Aid pouch per unit




+



X

Step 1: Choose Your Volume



Step 2: Choose Filtration

Hard to Filter (40 g DE/L)

Your Kit: **SDLV-0050-02E0-2**

1. 50 mL vacuum filtration unit
Sartolab® RF 50, 0.22 μm , PES
2. 2 x 1 g Filter Aid pouch per unit



+



X