




Application Recommendations

Reprocessing of Laboratory Glassware and Utensils

	neodisher® LaboClean FLA	neodisher® LaboClean FT	neodisher® LaboClean FM	neodisher® FA	neodisher® LaboClean A 8	neodisher® LaboClean F	neodisher® LaboClean LA	neodisher® LaboClean UW	neodisher® LaboClean GK	neodisher® LaboClean PLM	neodisher® MA	neodisher® EM	neodisher® Entschäumer S	neodisher® N	neodisher® Z	neodisher® LM 2	neodisher® LM 3	neodisher® LM 10	neodisher® PM 5
automated cleaning																			
mildly alkaline detergent									●	●	●								
alkaline detergent	●	●	●	●	●			●											
highly alkaline detergent						●	●												
universal applications	●				●														
cleaning enhancer												●							
defoaming agent													●						
acidic detergent														●	●				
neutralising agent														●	●				
manual cleaning																			
mildly alkaline detergent																●			●
alkaline detergent																	●	●	