

## For maximum freedom


- ✓ More work area in the laboratory: ideal for setting up analytical equipment with service connections
- ✓ Ceiling or wall mounted, depending on available space



# For your requirements

- ✓ Telescopic supports for flexible positioning
- ✓ Diverse service cassettes – suitable for your services
- ✓ Flexible cassette system enables custom fitting out and retrofitting





	Width			Boom height without suspension	Suspension length	Depth
Ceiling mounting	900	1200	1500	360	Adjustable up to max. 1214	200
Laboratory wall mounting	900	1200	1500	360		200

Dimensions [mm]