

IKA

designed for scientists



lab dancer

/// Data Sheet

Attractively designed, test tube shaker.

Designed for mixing small test samples with a touch function.

- Small, compact and reliable
- For small containers up to 30 mm in diameter, e.g. test tubes, centrifuge tubes, Eppendorf beakers
- Fixed speed (2800 rpm)
- Excellent mixing action
- The upper casing and the test tube surface are made from inert plastic. The bottom section is a coated, zinc die-casting

www.ika.com

Subject to technical changes



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide



designed for scientists

- Includes light 12 V power pack set





designed for scientists

Technical Data

| | |
|--|----------------|
| Type of movement | orbital |
| Shaker diameter [mm] | 4.5 |
| Permissible shaking weight (incl. attachment) [kg] | 0.1 |
| Motor rating input [W] | 1.2 |
| Motor rating output [W] | 0.8 |
| Permissible ON time [%] | 100 |
| Speed fix [rpm] | 2800 |
| Speed display | none |
| Timer display | none |
| Touch function | yes |
| Dimensions (W x H x D) [mm] | 100 x 70 x 100 |
| Weight [kg] | 0.776 |
| Permissible ambient temperature [°C] | 5 - 40 |
| Permissible relative humidity [%] | 80 |
| Protection class according to DIN EN 60529 | IP 40 |
| Voltage [V] | 100 - 240 |
| Frequency [Hz] | 50/60 |
| Power input [W] | 30 |

